

Automated sample preparation using a digital syringe with embedded SPE capability

Introduction

Popular sample preparation techniques such as solid phase extraction (SPE) can be time consuming and expensive, whether implemented on automated platforms or performed manually. Furthermore, method development to improve the efficiency of the sample preparation process is often difficult and time prohibitive, with long evaporation steps throughout the process.

A key objective for all laboratories is to implement efficient sample preparation where there is a real need to improve productivity and minimize waste. While the introduction of automated and semi-automated platforms into the laboratory for sample preparation enables improved efficiencies, they generally come at a significant cost.











We have combined the advantages of automation and SPE using a digital syringe with an embedded SPE cartridge. This not only improves efficiencies, but virtually eliminates solvent use and waste. The miniaturized format is ideal for small valuable samples such as biological extracellular fluids.

Here we demonstrate the advantages in combining the automation of a hand held, digital syringe with miniaturized SPE sorbent embedded in the needle of the syringe. Method development is rapid and inexpensive, enhancing laboratory workflow, while increasing the accuracy and reproducibility of the SPE process.

Background

MEPS® (Micro Extraction by Packed Sorbent) is a miniaturized version of SPE with several key advantages over traditional SPE (see Table 1). In its original format, MEPS is embedded in an analytical syringe and only suffers from the operator dependent steps of aspiration and dispensing where control of the flowrate typically varies from operator to operator.

Table 1. Comparison table SPE vs MEPS.

	SPE	MEPS
Sample volume	 3 mL	 50 µL
Time		
Cost		
Solvents	 10 mL	 500 µL
Evaporation step	 Yes	 Unnecessary

Coupling MEPS technology with a digitally controlled analytical syringe, eVol (Figure 1), enables semi-automated sample preparation which can be programmed to carry out the whole extraction through to injection into the analysis system.




Figure 1. Programmable digital analytical syringe - eVol XR, with eVol syringe, eVol MEPS syringe and MEPS BIN (Barrel Insert and Needle).

Method, results and discussion

Caffeine in saliva

Sample: Following caffeine consumption, saliva was collected into a clean glass vial. A 1 mL aliquot was diluted with an equal volume of saturated sodium tetraborate solution to buffer the sample to pH 9.5. The sample was then filtered and frozen for subsequent evaluation.

MEPS extraction: The extraction method was first optimized by investigating optimal flow rate, the number of extraction cycles and the number of elution volumes. A C18 MEPS BIN (Barrel Insert and Needle) and eVol MEPS syringe was coupled to the eVol, and used to extract caffeine.

The eVol program was divided into 8 functions (Condition MeOH, Condition H₂O, Load sample, Wash H₂O, pH adjust, Wash H₂O, Air dry, Elute) with a total of 41 steps which took around three minutes (see Table 2). The optimized method was repeated by two additional operators.


Table 2. eVol MEPS method for extraction of caffeine.

Step	Mode	Amount (µL)	Speed
1	Methanol		Prime
2	Aspirate	20	4
3	Dispense	20	4
4	Aspirate	20	4
5	Dispense	20	4
6	Aspirate	20	4
7	Dispense	20	4
8	H₂O		Prime
9	Aspirate	20	4
10	Dispense	20	4
11	Aspirate	20	4
12	Dispense	20	4
13	Aspirate	20	4
14	Dispense	20	4
15	Sample		Bind
16	Aspirate	50	4
17	Dispense	50	4
18	Aspirate	20	4
19	Dispense	20	4
20	Mix (x8)	50	4
21	H₂O		Wash
22	Aspirate	20	4
23	Dispense	20	4
24	Saturated sodium tetraborate		
25	Aspirate	20	4
26	Dispense	20	4
27	H₂O		Wash
28	Aspirate	20	4
29	Dispense	20	4
30	Air		Dry
31	Aspirate	50	4
32	Dispense	50	10
33	Aspirate	50	4
34	Dispense	50	10
35	Aspirate	50	4
36	Dispense	50	10
37	MeOH		Elute
38	Aspirate	20	4
39	Dispense	20	4
40	Aspirate	50	4
41	Dispense	50	10

Conditions

GCMS of the extracted caffeine from three operators performed on an Agilent 6890/573.



Column part number	054101
Phase	BPX5
Column	30 m x 0.25 mm x 0.25 µm
Oven	
Initial temperature	45°C for 4 min
Rate	10°C/min
Final temperature	300°C for 10 min
Injector	
Carrier gas	He
Carrier gas flow	1.2 mL/min
Constant flow	On
Injection mode	Splitless, 250°C
Purge flow	50 mL/min, nominal inlet pressure 127 kPa
Injection volume	1.0 µL
Liner type	FocusLiner with bottom taper
Detector	MS
MS transfer line temperature	280°C
Scan parameters	Low mass: 40.0, High mass: 500.0 Da at 2 scan/sec
MS quad temperature	150°C
MS source temperature	230°C



Opiate detection in urine




Normal urine contains a high concentration of salts and is not suitable to be injected straight into the MS. Ion suppression limits the detection of codeine and its metabolite morphine in the untreated sample.


A simple dilute and infuse experiment (Figure 3 and 4) demonstrates the inability to detect the ions of interest (476.19 Codeine-6-Glucuronide; 462.18 Morphine-6-Glucuronide; 286.15 Morphine and norcodeine; 272.13 normorphine).



eVol MEPS direct infusion into ESI

MEPS provides a fast and easy desalting step and the digital syringe format allows convenient sample introduction by direct infusion.





Conclusion

MEPS is a robust micro SPE technology suitable for high-throughput sample preparation applications.

The digital analytical syringe embodiment delivers controlled flow rates through the programmable drive, minimizing inter-operator error and enabling rapid method development.

The miniature format facilitates direct injection into the LC and/or GC, and as demonstrated here directly into the mass spectrometer source.

References

Esme Candish, Australian Centre for Research on Separation Science (ACROSS), School of Chemistry, University of Tasmania.

Information and support

Visit www.trajanscimed.com or contact techsupport@trajanscimed.com

Specifications are subject to change without notice.