

Accessories | Electronic battery

Electronic battery crimping tools effortlessly crimp/decap aluminum and 2-part aluminum/steel caps with just the push of a button.

- Consistently crimp caps to a reproducible seal setting with less physical strain.

Features and benefits

- Powerful brushless motor with a fast cycle.
- Long lasting lithium ion battery provides up to 800 crimps/decaps per charge with an average recharge time of 1-2 hours.
- Innovative charging system, which includes a warning light when recharging is necessary.
- Can be operated while charging.
- The controls to increase or decrease crimp settings are conveniently located on top of the unit, preventing accidental adjustment of the settings.
- OLED screen for simple operation displays tool size, battery status, crimp setting and cycle result.
- Performance tracking through an internal cycle log.

Recommended applications

High volume or strenuous crimping/decapping applications.

Product specifications

Tools in 8, 11, 13, and 20 mm crimping and decapping.



For more information about this product visit www.trajanscimed.com or contact techsupport@trajanscimed.com.