

## Accessories | Electronic high power unit

Ergonomic electronic high power crimping tools can effortlessly crimp/decap even the strongest all-steel caps with just the push of a button. The high power unit is designed as a single tool with interchangeable jawsets for the various functions and sizes needed in your laboratory.

- Consistently crimp caps to a reproducible seal setting with less physical strain.
- No need more multiple battery units - one high power unit covers all sizes and functions needed with the jawset attachments.

### Features and benefits

- Fixed power supply with cord (no battery).
- Uses multiple jawsets to achieve various functions and accommodate different size needs.
- The controls to increase or decrease crimp settings are conveniently located on top of the unit, preventing accidental adjustment of the settings.
- OLED screen for simple operation displays jawset size, crimp setting and cycle result.
- Performance tracking through an internal cycle log.
- Strong enough for all-steel caps.

### Recommended applications

High volume or strenuous crimping/decapping applications including all-steel caps and for requirements of multiple cap sizes and functions.

### Product specifications

Jawsets available in 8, 11, 13, and 20 mm crimping and decapping.

For more information about this product visit [www.trajanscimed.com](http://www.trajanscimed.com) or contact [techsupport@trajanscimed.com](mailto:techsupport@trajanscimed.com).

